

Oscillating Disk Rheometer



standard specification

Main power	220V, 50/60Hz, 5A Single phase
Oscillating frequency	100 cycle/minute
Oscillating angle	1', 3'
Temperature range	Natural ~ 200 °C
Temperature tolerance	±0.3 °C
Temperature control	Computer type, PID control
Temperature probe	PT 100Ω
Air pressure	Min. 4.2Kgf/cm ² (60psi)
Operating system	IBM P/C Compatible
Operating program for windows	
Data report	<ul style="list-style-type: none">●Max. Torque (MH)●Min. Torque(ML)●Initial cure time (ts1,ts2)●Cure time (tc10, tc50, tc90)
Statistics	Data files for EXCEL
Dimension	580 x 710 x 1168 (W x D x H)
Weight	230Kg

labkorea