

CENTRIFUGE for Personal LAB

CENTRIFUGE

Features

- Quiet, brushless DC motor
- Maintenance free



NB-501D

MICRO CENTRIFUGE



NB-501




NB-502

TABLE TOP CENTRIFUGE

SPECIFICATION				
MODEL	NB-501D	NB-501	NB-502	NB-504
SPEED RANGE	0~11000rpm	0~11000rpm	0~4000rpm	0~2000rpm
CONTROLLER	Micom, Key pad	Volume s/w	Micom, Key pad	Volume s/w
VORTEX				yes
DIMENSION	220x240x170mm	220x240x170mm	355x355x200mm	160x170x80mm
ROTOR(angle)	1.5mlx12	1.5mlx12	15mlx8	1.5mlx12
POWER	110/220V, 50/60Hz(select power)			

labkorea



NB-503



NB-503L

CONCENTRATOR

Features

- Quiet, brushless dc motor
- Maintenance free
- Precise vacuum sealing
- Unique coupling drive

SPECIFICATION		
MODEL	NB-503	NB-503L
SPEED RANGE	1000~4000rpm	1800rpm
DRIVE	B/L DC MOTOR	B/L DC MOTOR
CHAMBER	Ø150mm	Ø305mm
HEATER	ATC 60℃	ATC 60℃
VACUUM PUMP	Built in	
VACUUM GAGE	Built in	
MAX CAPACITY	1.5mlx12ea	1.5mlx60ea
DIMENSION	220x240x170mm	355x400x170mm
LID	Acrylic	Glass
POWER	110/220v, 50/60Hz(select power)	

labkorea