

# SHAKER for Personal LAB

## SHAKERS

### Features

Plate Type B/L Motor(microprocessor controlled speed)  
Constant speed under fluctuating power

- Orbital Shaker automatically regulates and stabilizes
- Flat, compact type shaker with a new oscillation mechanism
- Quiet running and precision circulation
- Variable attachments



NB-101S



NB-101SH



NB-101SS


### SPECIFICATION

MODEL	NB-101S	NB-101SS	NB-101SH
SPEED RANGE	25~350RPM		
STROKE	20mm		
PLATFORM SIZE	255x300mm(10"x12") with rubber		
OVERALL SIZE	280x330x70mm(11"x12.9"x2.75")		
HOLDER/RACK	Rubber pad	spring rack	9place x 250ml
WEIGHT	7Kg		
POWER	110/220V, 50/60Hz		

*labkorea*



NB-101MH



NB-101M




NB-4SM

**TABLE TOP SHAKER**

SPECIFICATION				
MODEL	NB-101M	NB-101MH	NB-101MS	NB-4SM
SPEED RANGE	25~350RPM			
STROKE	30mm(1.18")			
PLATFORM SIZE	460x455mm(18"x18") with rubber			4 STAGE
OVERALL SIZE	460x520x120mm(18"x20.5"x4.7")			500x560x180mm
HOLDER/RACK		23place 250ml	spring wire rack	
WEIGHT	35Kg	40Kg	40Kg	
POWER	110/220V, 50/60Hz(select power)			

OPTION ACC.



MOD. NO	MAX. CAPACITY	FLASK HOLDER
SH100	100ml	36ea
SH250	250ml	23ea
SH500	500ml	16ea
SH1000	1000ml	9ea
SH2000	2000ml	4ea
SSM	spring rack	

*labkorea*