

STIRRER for Personal LAB

STIRRER

Features

- Magnetic Stirrer
- Convenient use (immersion / dry)
- Separated speed controls
- Controlled with electric feedback(stabilization speed)




NB-401-I

IMMERSION STIRRER



NB-402S4

MULTI-POSITION STIRRER



NB-403S3

SPECIFICATION			
MODEL	NB-401-I	NB-402S4	NB-403S3
SPEED RANGE	35rpm~1500rpm		
STIRRING POSITION	1ea	4ea	3ea
PLATE MATERIAL	Stainless steel		
IMMERSION	yes		
DIMENSION	Ø140x60mm(h)	280x300x60mm	600x270x60mm
POWER	110/220V, 50/60Hz(select power)		

labkorea