

Multi Purpose Centrifuge



Model : MF 50 D

This is generally used for multi purpose and routine experiments related to urine, parasites, and blood. This reliable centrifuge has an outstanding price-performance ratio.

* Feature

- One chip microprocessor control
- Digital display of program, speed and time
- Equipped keyboard & alphanumeric display
- 10 program storage and recall
- Display and control RCF
- Gentle start-up and braking—three acceleration/deceleration profiles
- Interactive induction magnetic brake
- Set conditions can be changed while spinning
- Interlocking door for safety
- Display abnormality : imbalance, over-speed, door-open

* Specification

Max. RPM	6000 rpm
Max. RCF	4200 xg
Max. Capacity	15ml x 8
Temperature range	Air cooling
Timer	99min 59sec or free run
Braking rates	3 step
Programmability	10 memory
Drive Motor	High torque DC motor
Power	220Vac, 50/60Hz, Single phase, 1 Kw
Dimension	360W x 465D x 300Hmm
Net Weight	30 Kg

labkorea