


Homogenizer




- ❖ Homogenizing wastewater samples over the use in laboratory reactors.
- ❖ Mix emulsify your larger biological sample.
- ❖ Designed for tough solid or large volume samples disperse with stepless speed adjustment for quantities.
- ❖ Especially powerful disperse, due to high rotational speed (18,000 R.P.M)
- ❖ Also suitable for crushing of fibrous materials volumes.
- ❖ Easy exchange of dispersing element by quick disconnect coupling.
- ❖ Electronic safety circuit and smooth start.
- ❖ This unit includes motor, stand with clamp and tool kit.
- ❖ Digital tachometer separately main body (HMZ-20DN).
- ❖ Solve your problems with hard to wet fine powers and mixture solid.
- ❖ Very easy operation and anyone can use.
- ❖ Materials of solution-touch parts is stainless steel and Teflon.

TOOLS



- 3 : 2 Liter capacity
 - 2 : 5 Liter capacity
 - 1 : 10 Liter capacity
- with accessories



Specification

	HMZ-20N	HMZ-20DN	HMZ-C40	HMZ-30S
POWER	ONLY 110V			ONLY 220V
WATTAGE	200W			20W
R.P.M	18,000			12,000
R.P.M INDICATOR		DIGITAL DISPLAY	ALALOG DISPLAY	
VESSEL(STAINLESS STEEL)	x	x	300ML, 500ML, 1,000ML	x
SPEED CONT'	STEPLESS ELECTRONIC			
TIMER	OPTION		60min	

labkorea