

Hot Plate & Magnetic Stirrer



HMS-10D



HMS-10D



HMS-10



HMS-20



HMS-30

HMS-10,20,30

- Feature a ceramic-coated to plate for even heating.
- Pure ceramic top is acid and alkali resistance.
- The stepless heat adjustment let you dial in temperatures up to 370°C
- Quiet operation and regulation of motor speed by specially made DC motor.

SPECIFICATION

	HMS-10D	HMS-10	HMS-20	HMS-30
POWER	AC 220V			
MOTOR	DC MOTOR			
R.P.M	1,200 MAX	1,700 MAX	1,700 MAX	
SPEED CONTROL	FEEDBACK SPEED CONTROL	STEPLESS ELECTRONIC		
PLATE SIZE	180 x 180 (MM)			
TEMP CONTROL	P.I.D Type	PROPORTIONAL VOLTAGE CONTROL		PROPORTIONAL VOLTAGE CONTROL
TEMPERATURE	370°C MAX.			370°C MAX.
DIGITAL SETTING TEMP & SPEED	YES	NO	NO	NO

labkorea