


MIXER

SH2000



SH2000-DX







- Strong and stable velocity.
- Maximum flexibility

Rack changeable
STB20 - For mini orbital shaker

SPECIFICATION		
MODEL	SH2000	SH2000-DX
Orbit range	3mm	
Speed	Variable 0~2000rpm	
Timer	60min hold/dial	
Dimension(mm)	W170xD238xH178	W323xD286xH216
Temperature range	Ambient +5°C~95°C	
Temperature accuracy	±0.1°C	
Weight	9.8Kgs	14Kgs
Power	AC110V/45W(AC220V/45W)	AC110V/250W(AC220V/250W)

ACCESSORIES for sh2000

TUBE RACKS

MODEL	DEMONSTRATION	IMAGE
T0.5	0.5ml - 64holes	
T0.5,1.5	1.5ml - 32holes and 0.5ml - 32holes	
T1.5	1.5ml - 64holes	
STB20	Mini-Shaking plate	

labkorea