

Auto Ozone Tester



According to ASTM D 3395

A fully automatic Ozone Test Chamber, State of the art technology and construction techniques have been combined with expertise to produce a reliable, accurate, labour & cost saving unit, which is truly environmentally safe.

1. Standard Specification

- Main power : AC220V, 50/60Hz Single Phase
- Ozone range : Max.250pphm
- Ozone Tolerance : +/- 0.5%`C
- Temperature range : Max. 60`C
- Temperature Tolerance : +/-0.5`C
- Chamber size : 1,000 x 1,000 x 1,000(W x D x Hmm)
- Weight : 300kg

2. Control system

- Ozone monitor, controller and Generator
- Temperature controller
- Temperature sensor
- Recorder
- High Voltage Transformer

3. Dynamic System (Option)

- System : Motor, Cam and Sample Grip

4. Spare parts for two years operation

- 2 pieces of Ozone lamp
- 1 pieces of Heater
- 3 pieces of pen for recorder (blue, red color)
- 10 rolls of recording paper

labkorea