

DISINTEGRATION TESTERS

Sub disk mold

Reliable water bath accuracy

User friendly mechanism

Digital timer with LED display

Stroke height adjustable acc.to international pharmacopoeias



DIT-200



DIT-400

Specification

Model	DIT-200	DIT-400	DIT-600
Basket No.	2 ea	4 ea	6 ea
Motor	25W x 1	25W x 2	25W x 3
Stroke	30 cpm \pm 1 cpm 22,50,55,75mm		
Timer	Digital (99H 59M 59S) Elapsed time (with alarm)		
Temp.	30~45`C, 37`Csetting 0.1`C display & set		
Safety	ELBS/W & heater burn out protector		
Water bath	Stainless steel and safety glass		
Sub disc	12	24	36
Thermometer	3	5	7
Glass tube	12	24	36
Dimension (W x D x H)	400 x 295 x 450	700 x 295 x 450	1000 x 295 x 450
Capacity	10 L	18 L	27 L
Power	AC 110,220V 50/60Hz		
Weight	15 kg	30 kg	40 kg

* Designed acc.to USP, BP, JP & other Pharmacopoeias

* Sub basket for Granule is accessories

MODEL

DIT - 200

DIT - 400

DIT - 600

labkorea