

DISSOLUTION TESTER

Best Solutions

* Main body & water bath



* Accessories



Microprocessor control and digital display

Reliable servo motor

Ergonomic vessel design

Automatic lifter

Time monitor

Specification

Model	DST-610	DST-810
Number of shaft	6 ea	8 ea
Speed range	2 - 200 rpm	
Accuracy of speed	± 1 rpm	
Drive system	Servo motor	
Temp. control	Digital auto control 0.1 ± 1 °C	
Water bath	Stainless steel & safety glass	
Circulation pump	1/20 Hp	
Time monitor	Digital display time & alarm	
Speed monitor	Digital display (0-250rpm)	
Dimension(overall)	600W x 465D x 650Hmm	
Power	AC 110/220, 50/60Hz	
Vertical lifter	Automatic Motor driven	
Rotating basket	6 ea	8 ea
Paddle	6 ea	8 ea
Sinker	6 ea	8 ea
Vessel (with lid)	7 ea	8 ea
Power consumption	1KW	1.2KW

* All models designed acc.to USP, JP & other Pharmacopoeia...

MODEL DST-610 DST-810

labkorea