


FRIABILATORS

Reliable geared motor & synchronous belt

Digital timer & alarm buzzer

Friability test drum

Abrasion test drum



MODEL FAT-10 FAT-20

Specification

Model	FAT - 10	FAT - 20
Number of drum	1 ea	2 ea
Motor	Synchronous motor 15W	
RPM	25 rpm	
Timer	Digital (0~59min. 59sec.)	
Drum size	300mm dia. X 35mm width clear acryl	
Dimension (overall)	300W x 250D x 340mm	420W x 250D x 340Hmm
Power consumption	AC100,220V 50/60Hz	
Weight	8 kg	10kg

labkorea