

ROTATOR



AG

- Strong and stable velocity
- Flexible angle of rotation plate (from 0~90)

SPECIFICATION	
MODEL	AG
Motor	.DC-3.5A .Engineering plastic gear
Angle range	0°~ 90°
Speed	Variable 0~50rpm
Timer	3hrs with hold / dial
Temperature resist	Lower then 65°(dry air)
Dimension(mm)	W220×D270×H245 (W8.7"×D10.6"×H9.6")
Weight	7.5kgs (16.5 Lbs)
Power	AC 110V (AC 220V)

ACCESSORIES For AG (ROTATOR)

D1.5-32	1.5ml Eppendorf Tube	32ea	
D1.5-60	1.5ml Eppendorf Tube	60ea	
D15-16	15ml conical Tube	16ea	
D50-8	50ml conical Tube	8ea	
D250-4	250ml Flask	4ea	

MODEL	SIZE	COMENT
BT40-120	40mm diameter X 120mm length	with silicon cap
BT40-200	40mm diameter X 200mm length	with silicon cap
BT60-200	60mm diameter X 200mm length	with silicon cap
BT80-200	80mm diameter X 200mm length	with silicon cap






HAG






- Strong and stable velocity
- Flexible angle of rotation plate (from 0~90)
- Rocking platform and rotor could be attached.

SPECIFICATION	
MODEL	HAG
Motor	.DC-3.5A .Engineering plastic gear
Angle range	0°~ 90°
Speed	Variable 0~50rpm
Timer	3hrs with hold / dial
Temperature resist	Lower then 65°(dry air)
Dimension(mm)	W470×D270×H245 (W18.5"×D10.6"×H9.6")
Weight	8kgs (17.6 Lbs)
Power	AC 110V (AC 220V)


ACCESSORIES For HAG (ROTATOR)

MODEL	BOTTLE SIZE	CAPACITY	IMAGE
RT-4	80mm diameter bottle	4ea	
RT-6	60mm diameter bottle	6ea	
RT-10	40mm diameter bottle	10ea	

RT-8V	50mm conical Tube	8ea	
RT-24H	50mm conical Tube	24ea	
RT200-64	1.5ml Eppendorf Tube	64ea	
RT200-96	1.5ml Eppendorf Tube	96ea	
RT200-128	1.5ml Eppendorf Tube	128ea	
RT300-120	1.5ml Eppendorf Tube	120ea	
RT300-180	1.5ml Eppendorf Tube	180ea	
RT300-240	1.5ml Eppendorf Tube	240ea	

D1.5-32	1.5ml Eppendorf Tube	32ea	
D1.5-60	1.5ml Eppendorf Tube	60ea	
D15-16	15ml conical Tube	16ea	
D50-8	50ml conical Tube	8ea	
D250-4	250ml Flask	4ea	

RK-H



MODEL	SIZE	COMENT
BT40-120	40mm diameter X 120mm length	with silicon cap
BT40-200	40mm diameter X 200mm length	with silicon cap
BT60-200	60mm diameter X 200mm length	with silicon cap
BT80-200	80mm diameter X 200mm length	with silicon cap

labkorea