

# THERMO REGULATOR



## TR500

- Ordinary container will be changed into an accurate water bath by this thermo-regulator.
- Thermo-controller has following functions auto-tuning, solid state-proportional control and limitation of setting temperature.
- As the thermo-controller is easy to assemble and very convenient to maintain.

### SPECIFICATION

MODEL	TR 500	
Temperature controller(digital)	temp. control range	ambient.+5°C~100 °C(additionalath:70°C )
	temp. uniformity	better than $\pm 0.1^{\circ}\text{C}$
Clamping thick.	25mm(0.95")(bath well)	
Dimension(mm)	out side	W124× D175× H287(W4.9"× D6.9"× H1.3")
	water bath in side	W124× D140× H127(W4.9"× D5.5"× H5")
Weight	3.5kg(7.7Lbs)	
Power	AC 110V	

*labkorea*