

# ULTRASONIC

C L E A N E R



Powerful and reliable cleaning  
Versatile performance  
- heating, timer  
Ergonomic design  
Durable structure



**A** powerful industrial BLT trembler is adapted to produce a high efficient oscillation. It ensures powerful and reliable cleaning. The output can be set to high, medium, low in 3 stages. You can change the ultrasonic power for your application.

**T**he silicone rubber heater makes the cleaning liquid heat up to the temperature required, so that the efficiency can be increased. Also, you can set the timer up to 99 minutes and the remaining time is displayed. Control panel and temperature and time displays are designed ergonomically.

**E**xternal case and cover made of waterproof and corrosion resistant plastic are rigid and durable. A Basket is included as a standard and can be comfortable to handle the sample.



## 4020, 2010, 1505 Ultrasonic Cleaner

### Specification

Model	4020(400W) 4020(500W)	2010(200W) 2010(300W)	1505(150W) 1505(200W)	1002(100W)
Capacity	20 L	10 L	5.7 L	2L
Temp. range	0-70 ; É			
Timer	99 min.			
Transducer	BLT			
Ultrasonic frequency	40KHz			
Ultrasonic power	400W 500W	200W 300W	150W 200W	100W
Body material	ABS			
Tank Size(WxDxH)	450x290x160 mm	295x240x160 mm	295x150x160 mm	142x132x104mm
External size(WxDxH)	721x451x297 mm	490x368x297 mm	455x298x267 mm	*
Electric requirement	AC 110V, 50/60Hz or AC220V, 50/60Hz			
Weight(Net)	17.0 Kg	11.0 Kg	8.0 Kg	*

*labkorea*